


MUSE™ System - Transoral Fundoplication for GERD



 **MUSE™** SYSTEM
GERD Relief in Sight

Medigus Ultrasonic Surgical Endostapler

Distal Tip



Medigus Ltd. is a leading international company specializing in minimally invasive medical procedures.

The revolutionary MUSE™ System enables a single endoscopist to perform a transoral anterior fundoplication for the treatment of GERD, avoiding abdominal incisions.

System components:

- A flexible surgical endostapler:
 - Standard titanium surgical staples
 - Video capabilities
 - Ultrasonic sight
- Digital console

Key features:

- Transoral procedure - No scars
- Video and ultrasound guided surgical endostapler - Fully controlled procedure
- Standard surgical stapler - Reliable fundoplication
- Intuitive operation by a single endoscopist - Fewer resources needed

The MUSE™ System has FDA clearance for marketing and CE mark.


- ✦ Video and ultrasound guided surgical endostapler
- ✦ Standard surgical staples
- ✦ Intuitive operation by a single endoscopist



MUSE™ System Procedure – Step by Step


Advance into stomach and retroflex

1



Retract MUSE™ System to 3cm proximal to GEJ, clamp tissue and staple

2



Repeat to create flap-valve

3

